



MOH-ZNPH DIGITAL POE SCREENING MOBILE APP SCREENER'S MANUAL

INSTALLATION




How to Install the App on Your Android Device . Please follow these step-by-step instructions to complete the installation.

- **Step 1:** Download the App from **WhatsApp/Provided link and Open it.** You might see a notification indicating the file type. If any warnings appear, you can safely ignore them and tap “OK” to proceed. If prompted to allow your phone to scan the file, select Scan. If permission requests are shown, follow the on-screen instructions to grant the necessary permissions.
- **Step 2:** Enable Installation from Unknown Sources. Since the app is not available on the Google Play Store, you must enable installations from unknown sources. Follow these internationally recommended guidelines: Access Device Settings: After downloading the APK, you may receive a notification that the installation is blocked due to security settings.
- **Enable Unknown Sources:** Open the Settings app on your Android device. Navigate to Security. (On some devices, you might need to go to Apps & Notifications > Special App Access > Install Unknown Apps.) Locate your messaging app (e.g., WhatsApp) and enable Allow from this source. Note for Android 8.0 (Oreo) and later: Installation permissions are managed on a per-app basis, so ensure that the permission is enabled specifically for the app you’re using to install the APK.
- **Install the APK:** Once the setting is enabled, tap the downloaded Digital POE Screening.apk file in your notifications or file manager to start the installation. Follow the on-screen instructions to complete the process.

INSTALLATION

- Step 3: Install the App Locate the APK File: Go to your Downloads folder using a file manager or tap the notification when the download completes. Start Installation: Tap on the APK file. You may be asked to confirm if you really want to install the app; just click Install. Wait for the Installation: The app will install in a few seconds.
- Step 4: Open the App Once the installation is complete, you can click Open to start using the app. The app icon will also be added to your home screen or app drawer for easy access later.

INSTALLATION



POE Health Screening

Digital Tool

Username

Password

LOGIN

LOGIN

Enter Username and Password

Username : ecsa

Password 12345678




POE Health Screening

Digital Tool

Username

Password

LOGIN



LANDING PAGE AFTER LOGIN IS THE PRIMARY SCREENING DASHBOARD

- This shows primary screening analytics for head count and gender
- It allow filtering of the analytics by date , month and year

Search menu items

Rwanda POE Screening

Welcome, ecsa

Primary Screening

Secondary Screening

Dashboard

Management

POE Management

User Management

Secondary Screening

Secondary Dashboard

Secondary Screening

Primary Dashboard

Session Management

Logout

PRIMARY SCREENING

Screening for people who do not present any concerning symptoms or signs or have not reported any health complaints. It records gender and head count

PRIMARY SCREENING FORM

The image displays two mobile application screens for a 'Primary Screening' form.


Left Screen (Initial Form):

- Header:** 'Primary Screening' with a menu icon and a refresh icon.
- Title Bar:** 'Primary Screening Form' with a refresh icon.
- Question:** 'Does the traveler have any of the following symptoms: fever (38°C or higher), feeling unwell, continuous cough, difficulty breathing, diarrhea, vomiting, skin rash, chills, unusual bruising or bleeding, limb weakness, dizziness, or jaundice?' (This text is on a blue background in the image).
- Data Sync Status:** A checkmark icon and the text 'Data Sync Status Synchronized with Central Database'.
- Symptoms:** A toggle switch for 'Does the traveller have any symptoms?'.
- Gender:** A dropdown menu currently showing 'Male'.
- Buttons:** A blue button labeled 'No Symptoms, Proceed' with a person icon, and a 'Previous' button with a left arrow.

Right Screen (Screening Result):

- Header:** 'Primary Screening' with a menu icon and a refresh icon.
- Title Bar:** 'Primary Screening Form' with a refresh icon.
- Section Header:** 'Screening Result'.
- Classification:** A green box labeled 'Classification Non-Case'.
- Recommended Action:** A yellow box labeled 'Recommended Action Issue safety guidelines and clear for travel.'.
- Buttons:** A green button labeled 'Complete Screening' with a checkmark icon, and a 'Previous' button with a left arrow.

- If a traveler does not have any of the following : fever (38°C or higher), feeling unwell, continuous cough, difficulty breathing, diarrhea, vomiting, skin rash, chills, unusual bruising or bleeding, limb weakness, dizziness, or jaundice? , The screener continues to record gender and complete the primary screening
- For primary screening to conclude successfully, The green complete screening button at the end of the form has to be clicked



SECONDARY SCREENING

Digitized version of the national POE case investigation form

HOW TO TRIGGER SECONDARY SCREENING

- Once a traveler has any of the indicated signs or symptoms
- Symptoms button is turned on and a form appears, Enter the traveler name who is being referred to secondary screening and click proceed
- Finally click the refer to secondary screening button and a notification for a traveler pending secondary screening will be created. This traveler is then moved to the holding area to continue with secondary screening

Traveller Information

Close

Traveller Name

Proceed

Primary Screening

Primary Screening Form


Screening Result

Classification
Suspected Case

Recommended Action
Refer to secondary screening for further examination.

Refer to Secondary Screening ↗

← Previous



START SECONDARY SCREENING

- On the app menu , Select secondary screening , You will now be taken to all travelers pending secondary screening

Secondary Screening

Referral Alerts

Q Search by name or ID

ahah

ID: SCR-921155

Created: 10/14/2024, 2:51:47 PM

pending

Christopher Minja

ID: SCR-186338

Created: 10/14/2024, 2:51:47 PM

pending

gato boy

ID: SCR-245232

Created: 10/14/2024, 2:51:47 PM

pending

gato boy

ID: SCR-825425

Created: 10/14/2024, 2:51:47 PM

pending

gato

ID: SCR-336092

Created: 10/14/2024, 2:51:47 PM

pending

gato

ID: SCR-812766

Created: 10/14/2024, 2:51:47 PM

pending

Secondary Screening

Referral Details

CLOSE

ahah

ID: SCR-921155

Symptoms

Traveller reported symptoms

Status

pending

Referred By

ecsa (admin)

ict@ecsahc.org

Point of Entry

Rusumo Border Post (land_border)

Kirehe, Eastern Province, Rwanda

Created At

10/14/2024, 2:51:47 PM


START SECONDARY SCREENING

CANCEL REFERRAL

SELECTING A TRAVELER

- From the secondary screening notification select a traveler to screen
- Upon selection , The referral details are shown in a pop up,
- Click start secondary screening to start screening a patient or cancel referral to delete the notification if applicable

12

Secondary Screening 

Open Secondary Screening Form

Welcome to the Rwanda POE secondary screening process. Please begin by syncing data.


Start Screening →

SECONDARY SCREENING FORM

1. Now fill in the whole form by clicking start screening
(Symptom assessment, Traveler information step 1- 6)
2. Countries should be selected from the drop-downs NOT manually typed for fields that require countries


COMPLETE SECONDARY SCREENING

- Upon completion of the form
- Case classification and recommended course of action is presented to the screener
- All cases that require referrals should have their referral data recorded first by clicking refer to nearest health facility before completing screening
- App case recommendation and course of action should be read and understood by the screener
- Click complete screening to create a national alert and record case




POE MANAGEMENT


1. We need to capture the referral form and complete it.



USER MANAGEMENT

- Create, Update and Delete users/screeners/admins
- Begin by selecting the POE whose users are to be managed





SECONDARY SCREENING DASHBOARD

- Shows analytics about secondary screenings done at the POE
- Reports are filtered by time
- Screened cases can also be seen in detail by selecting a case from recent screenings

An abstract geometric design on the left side of the slide. It features a dark blue background with various geometric shapes and patterns. A white circle is positioned near the top left. Below it, a light blue semi-circle is visible. To the right of the semi-circle, there is a pink triangle with diagonal lines. Below the semi-circle, there is a pink square with a pattern of concentric lines. To the right of the square, there is a light blue triangle. Below the square, there is a pink triangle. To the right of the triangle, there is a dark blue triangle. The overall design is modern and minimalist.

THANK YOU